



# Activ+ Mattress

## Comfort Range

Experience everyday comfort and support with the Activ+ Mattress Comfort Range, thoughtfully designed to adapt to your body, keep you cool, and move effortlessly with adjustable beds.



### Designed for Adjustable Beds

Dual-notched foam layers and sidewall hinges allow the mattress to flex seamlessly with adjustable bed functions while maintaining comfort and support

### Waveform Core

Designed with varying wave frequencies to provide targeted comfort and support where it is needed most. Firmer support zones at the head and foot sections help promote stability, while the torso zone encourages natural body alignment

### Australian-made Quality


Proudly made in Australia using medical-grade, high-density foam for lasting comfort, support, and durability

### Open Cell Foam

A highly breathable surface layer promotes airflow throughout the mattress, helping to release excess heat and maintain sleeping comfort

### Hybrid Foam

Adapts to body contours for personalised comfort and support, while its responsive design assists with repositioning and in-bed mobility



### Anti-bacterial, Stretch-knit Covers

Standard Knit, Cooling Gel, and Pressure Care

30mm Breathable Top Hybrid Foam Layer

Reinforced Side Layer Hybrid Foam

60mm Mid Layer Notched Hybrid Foam

60mm Base Layer Notched Hybrid Foam

### Anti-bacterial, Stretch-knit Covers

The mattress covers are made with silver-infused fibres, commonly used to help maintain a fresh and hygienic sleep surface. The range is available with either a Stretch-Knit Cover, offered with Classic and Cooling Gel interchangeable cover tops, or a Pressure Care cover designed for pressure care applications.

### Stretch Knit Covers *Interchangeable Covers for year-round comfort*



#### Classic

Balanced support for all users



#### Cooling Gel

Promotes a cool, dry sleeping climate

#### Soft

SWL 175kg

BEM469950 - Long Single  
BEM469956 - King Single  
BEM469962 - Long Double  
BEM469968 - Queen

#### Medium

SWL 200kg

BEM469952 - Long Single  
BEM469958 - King Single  
BEM469964 - Long Double  
BEM469970 - Queen

#### Firm

SWL 250kg

BEM469954 - Long Single  
BEM469960 - King Single  
BEM469966 - Long Double  
BEM469972 - Queen

### Pressure Care Covers

#### PU Cover

Designed for moisture management and protection in pressure care environments

#### Soft

SWL 175kg

BEM469974 - Long Single  
BEM469980 - King Single  
BEM469986 - Long Double  
BEM469992 - Queen

#### Medium

SWL 200kg

BEM469976 - Long Single  
BEM469982 - King Single  
BEM469988 - Long Double  
BEM469994 - Queen

#### Firm

SWL 250kg

BEM469978 - Long Single  
BEM469984 - King Single  
BEM469990 - Long Double  
BEM469996 - Queen

**7 YEARS LIMITED WARRANTY\***

\*Foam 7 years | Stretch Knit Covers 2 years | PU Pressure Reducing Cover 5 years