



AEROLITE CARBON FIBRE POWERCHAIR

The AeroLite Powerchair is light... VERY light. The modern carbon fibre frame is so light that it is a cinch to fold, lift and carry, to easily load into a car or to store when not in use.

AeroLite's lithium-ion batteries are approved for air travel and yes it comes with two of them! Perfect for taking you further for longer outings with no airport stress and less hassle and worry about battery charging down-time. Intuitive joystick control, contoured comfort seating and puncture-proof tyres make for confident travel on different surfaces, indoors and out.

The AeroLite makes everyday life easier for users and carers alike.





Foldable Design
Ideal for travel and easy
to lift into the car

13.2kg
Minimum
Frame Weight

20km
Range

Two
batteries
included

SWL
120kg
19st

2 YEAR WARRANTY
Frame

1 YEAR WARRANTY
Components

Australia
1300 133 120 | aspirecare.com.au

United Kingdom
01622 541235 | aidacare.co.uk

Key information about the Aerolite Carbon Powerchair



The intuitive joystick simplifies navigation, reducing strain and promoting greater comfort and efficiency in daily activities.


Accessories




Cup Holder
SCA669200




Phone Holder
SCA669210




Tool-free assembly and disassembly make setup quick and hassle-free



Two lithium-ion batteries are included providing a leading 20km driving range (1 battery to operate and 1 spare)
IATA-approved for air travel**



Low-maintenance, puncture-proof solid tyres let users travel confidently on different surfaces



Flip-up footplates improve transfer safety during seat entry and exit

CODE: PWS654150

Specifications	mm	in
Overall L x W x H	960 x 630 x 895	38 x 25 x 35
Folded L x W x H	960 x 315 x 800	38 x 12 x 31
Seat Depth	400	16
Seat Width	450	18
Anterior Seat Height	480	19
Posterior Seat Height	440	17
Backrest Height	420	17
Backrest Width	425	17
Armrest Height	235	9
Armrest Length	350	14
Max Obstacle Height	25	1
Min Turning Radius	900	35
Max Climbing Angle	6°	
Max Driving Speed	6km/h	
Driving Range	20km	
Front Wheels	180	7
Rear Wheels	300	12
Electronics		
Battery (Lithium-ion)	1 x 24V 12Ah (plus 1 x 24V 12Ah spare)	
Battery Charger	2.6A	
Motor/Drive System	150W	
Weights	kg	lbs
Min Product Weight (No Battery or Foot Plate)	13.2	29
Total Weight	16.6	36.7
Battery Weight	1.45kg (each)	3.2
Footplate Weight	0.5	1.2

*Range of motion dimensions are accurate to a +/-10mm tolerance
**The AeroLite's Lithium-Ion battery is compliant with international regulations. Please check with airlines for specific compliance regarding UN3171.