



# DELTA SCOOTER

## SMALL PULL APART

**Smart Design. Strong Mobility. Travel Ready.**

Engineered for convenience and performance, the Aspire Delta Small Scooter combines durable construction with practical features that make everyday mobility easier and more enjoyable both at home and on the go.





### Class-leading 150kg SWL

Offers enhanced support and stability for a wider range of users.



### Adjustable Tiller for Effortless Control

Easily adjusts to your preferred angle for a more comfortable and personalised ride.



### Compact Pull-apart Design

Easily disassembles into five parts for effortless lifting and transport.



### Additional Battery Options Available

Standard range battery up to 24km. Extended and airline-approved batteries also available (up to 35 km).



### 4-wheel Independent Suspension

Enhanced suspension absorbs bumps and uneven surfaces, for a smoother ride.



### Tight Turning Radius

Experience seamless navigation in tight spaces, perfect for confident indoor mobility.



### Smooth Transfers, Ultimate Comfort

Featuring a height-adjustable, swivelling padded seat, this scooter ensures easy side transfers every time.



### Intuitive Controls

Navigate with confidence using an intuitive dashboard featuring an LED battery display, speed controls, and a horn.

|                                    |                        |
|------------------------------------|------------------------|
| Overall L x W x H                  | 1048 x 550 x 900mm     |
| Seat Size D x W                    | 400 x 450 - 550mm      |
| Seat to Footrest Height            | 400 - 475mm            |
| Armrest Height                     | 180mm                  |
| Min Turning Radius                 | 1200mm                 |
| Ground Clearance                   | 70mm                   |
| Max Climbing Angle                 | 6°                     |
| Max Driving Speed                  | 6.4km/h                |
| Driving Range                      | 24km (35km optional)   |
| Front/Rear Wheels (puncture proof) | 190mm                  |
| Battery                            | 1 x 24V 20Ah (Lithium) |
| Product Weight (excl. battery)     | 35.6kg                 |
| Heaviest Part                      | 15.7kg                 |
| Battery Weight                     | 3.2kg                  |
| Motor                              | 270W, 2 Pole, Brushed  |
| Controller                         | PG 45A                 |
| Brake System                       | Electromagnetic        |
| Free Wheel Mode                    | Yes                    |
| Suspension                         | Yes                    |
| Charger Size                       | 24V 2A                 |

\*Scooter Frame/Motors/Axles/Brakes 2 Years | Battery 6 Months | Electronics 1 Year



SCT671600  
**Red**



SCT671610  
**Champagne**

### Optional Batteries

SCA669400  
**Airline-approved battery (11.2km) 24V 12Ah**

SCA669410  
**Extended range battery (35km) 24V 30Ah**

**2** YEAR WARRANTY\*

