



# ELIO SCOOTER

## SMALL PULL APART

### Smart, Strong and Exceptionally Portable

Experience the perfect blend of strength, comfort and convenience with the Aspire Elio Small Scooter. Built for an active lifestyle, Elio's lightweight design, flexible battery options, and comfort features offer confident, independent mobility, perfect for shopping, social visits, or days out.





### Light to Lift

At just 10.8kg, the heaviest part is easy for carers to lift and fits neatly in most car boots or storage spaces.



### Ultra-lightweight Aluminium Frame

Easily disassembles into five compact parts for effortless transport and storage.



### Additional Battery Options Available

Standard range battery up to 24km. Extended and airline-approved batteries also available (up to 35 km).



### Small Turning Circle

Ideal for navigating through tight spaces such as hallways, lifts, or small shops without frequent repositioning or reversing.



### Shock-absorbing Suspension

4-wheel independent suspension absorbs bumps from uneven surfaces for a smooth ride.



### Intuitive Controls

Navigate with confidence using an intuitive dashboard featuring an LED battery display, speed controls, and a horn.



### Smooth Transfers, Ultimate Comfort

Featuring a height-adjustable, swivelling padded seat, this scooter ensures easy side transfers every time.



### Flexible Tiller Angle for Personalised Handling

Quickly adjusts to your ideal position, delivering a smooth and comfortable driving experience.

|                                    |                        |
|------------------------------------|------------------------|
| Overall L x W x H                  | 1048 x 550 x 945mm     |
| Seat Size D x W                    | 400 x 450 - 550mm      |
| Seat to Footrest Height            | 400 - 475mm            |
| Armrest Height                     | 180mm                  |
| Min Turning Radius                 | 1200mm                 |
| Ground Clearance                   | 60mm                   |
| Max Climbing Angle                 | 6°                     |
| Max Driving Speed                  | 6.4km/h                |
| Driving Range                      | 24km (35km optional)   |
| Front/Rear Wheels (puncture proof) | 190mm                  |
| Battery                            | 1 x 24V 20Ah (Lithium) |
| Product Weight (excl. battery)     | 28.6kg                 |
| Heaviest Part                      | 10.8kg                 |
| Battery Weight                     | 3.2kg                  |
| Motor                              | 270W, 2 Pole, Brushed  |
| Controller                         | PG 45A                 |
| Brake System                       | Electromagnetic        |
| Free Wheel Mode                    | Yes                    |
| Suspension                         | Yes                    |
| Charger Size                       | 24V 2A                 |

\*Scooter Frame/Motors/Axles/Brakes 2 Years | Battery 6 Months | Electronics 1 Year



SCT671700  
**Blue**



SCT671710  
**White**

### Optional Batteries

SCA669400  
**Airline-approved battery (11.2km) 24V 12Ah**

SCA669410  
**Extended range battery (35km) 24V 30Ah**

**2** YEAR WARRANTY\*

