




EVOKE 2 WHEELCHAIR

ATTENDANT BRAKE OPTIONS

Convert your Evoke 2 wheelchair into a transit wheelchair with 12in wheels and attendant brakes or update your Evoke 2 with 24in rear wheels to include attendant brakes with these two accessory kits.



Evoke 2 Attendant Brakes Kit

WCA848592 - No Wheels

WCA848594 - Brakes
and Wheels



UPHOLSTERY
Silver/Black

FRAME
Silver



1 YEAR WARRANTY
Components

2 YEAR WARRANTY
Frame



ALUMINIUM	Seat Width (mm)	Overall Width (mm)	Seat Depth (mm)	Seat Height (mm)	Backrest Height (mm)	Armrest Height (mm)	Min. Chair Weight (kg)
MWS449500	350	560	300-450	440-480	430-490	205-225	14.8
MWS449510	400	610	300-450	440-480	430-490	205-225	14.9
MWS449520	450	660	300-450	440-480	430-490	205-225	15.0
MWS449530	500	710	400-550	440-480	430-490	205-225	15.1