



EVOKE 2 WHEELCHAIR

The Aspire EVOKE 2 and EVOKE 2 HD are highly adjustable, lightweight, aluminium wheelchairs designed to enable a personalised fit for every user. A comprehensive range of compatible accessories is available for increased independence and enhanced clinical outcomes.



CRASH TEST
APPROVED
AS/NZS 3696.19

1 YEAR WARRANTY
Components

2 YEAR WARRANTY
Frame



SILVER
FRAME



BLACK
UPHOLSTERY

EVOKE 2

SWL
150kg

EVOKE 2 HD

SWL
190kg

1300 133 120 | aspirecare.com.au

Swing-In / Out, Height Adjustable and Removable Foot Rests: This facilitates user transfers and allows optimal positioning for the changing needs of the user


Strong, Nylon Upholstery with a Comfortable Cushion and Padded Backrest: For increased comfort, a padded overlay is available as an accessory

Depth Adjustable Seating with Rear Angle Adjustability: For more accurate positioning, customisability and increased user comfort

Quick Release 24 inch Solid Rear Tyres: The wheels are easily removed to reduce the weight of the chair for transport and storage

Pictured with optional transit conversion kit

ADJUSTMENTS



View our Evoke 2
Wheelchair Attendant
Brake Options



View our
Aspire Wheelchair
Accessories Range



EVOKE 2 WHEELCHAIR	Seat Width mm	Seat Depth mm	Seat Height mm	Overall Width mm	Backrest Height mm	Armrest Height mm	Min Chair Weight kg	SWL kg
MWS449500 Standard 350mm	350	375 – 450	450 – 500	560	430 – 490	210 – 230	14.8	150
MWS449510 Standard 400mm	400	375 – 450	450 – 500	610	430 – 490	210 – 230	14.9	150
MWS449520 Standard 450mm	450	375 – 450	450 – 500	660	430 – 490	210 – 230	15	150
MWS449530 Standard 500mm	500	400 – 550	450 – 500	710	430 – 490	210 – 230	15.1	150
MWS449540 HD 570mm	570	400 – 550	450 – 500	780	430 – 490	210 – 230	16	190
MWS449550 HD 610mm	610	400 – 550	450 – 500	820	430 – 490	210 – 230	16.2	190

Note: Range of motion dimensions are accurate to a +/-10mm tolerance