




# HYBRID PATIENT SLING



## SLING SIZING

Coloured binding tape identifies sling size



**1** YEAR WARRANTY



MESH FABRIC



SMALL & MEDIUM



LARGE & X LARGE

- **Hybrid sling** facilitates both hygiene and general purpose manual handling
- **Breathable mesh** promotes air circulation, reducing moisture and heat build-up while helping to maintain skin integrity
- **Intuitive shape** allows ease of access for toileting and dressing
- **Security belt and padded leg straps** support the trunk and thighs
- **Full coverage** shoulder and head support
- **Durable loop webbing** designed not to tangle or 'cork-screw' during cleaning
- **Embroidered identification labelling** ensures equipment traceability, supporting inventory control and maintenance protocols
- **Machine washable** up to 85°C

CODE	DESCRIPTION	EDGING COLOUR	SWL
LSS390850	Aspire Hybrid Patient Sling Loops - Small	●	200kg
LSS390852	Aspire Hybrid Patient Sling Loops - Medium	●	200kg
LSS390854	Aspire Hybrid Patient Sling Loops - Large	●	300kg
LSS390856	Aspire Hybrid Patient Sling Loops - X Large	●	300kg