



JUNO SCOOTER

LARGE

Smooth, Confident, and Comfortable. Wherever Life Takes You.

Rediscover the joy of getting out and about with a scooter designed for all-day comfort and effortless rides. From smooth suspension to a plush captain's seat, every detail is built to make your journey more enjoyable, whether you are on local paths, shopping trips, or uneven outdoor terrain.





Reliable Performance

High-torque motor (850W) allows users to navigate steeper inclines.



Plush, Supportive Seating

Features a deeply padded, adjustable captain's seat with suspension, and flexible tiller. Seat slides, swivels and reclines for perfect positioning.



Smart Safety

Emergency stop button and automatic speed reduction while turning keeps you safe and secure.



Comfortable Ride

Independent suspension absorbs bumps and reduces vibrations for comfort on uneven surfaces.



Digital Dashboard

Clearly shows battery life, speed, and other key information, giving you confidence while you travel.



Effortless Steering in Compact Spaces

A 1660mm turning circle allows smooth manoeuvring around shop aisles and community spaces.

Overall L x W x H	1490 x 720 x 1400mm
Seat Size D x W	440 x 530-580mm
Seat to Deck Height	470mm
Seat to Ground Height	680mm
Armrest Height	180-270mm
Min Turning Radius	1660mm
Ground Clearance	100mm
Max Climbing Angle	10°
Max Driving Speed	10km/h
Driving Range	40km
Tyres	330mm Pneumatic
Battery	2 x 12V 75Ah (Lead Acid)
Product Weight (excl. battery)	96kg
Battery Weight	51.5kg
Power Motor	850W 4 Pole
Controller	RHINO 2 160A
Brake System	Electro-Mechanical
Free Wheel Mode	Yes
Charger Size	24V 8A



SCT671296
Storm Grey



SCT671295
Desert Sand



*Scooter Frame/Motors/Axles/Brakes 2 Years | Battery 6 Months | Electronics 1 Year

2 YEAR WARRANTY*

SWL
181kg

ASFY185.5.170925