

# DA VINCI QUATTRO ACTION

## 4 MOTOR LIFT RECLINE CHAIR

This modern lift recline chair combines quattro-action positioning with integrated lift and recline assistance, featuring pressure-redistributing chaise design and ergonomic lumbar support for optimal user comfort and enhanced sit-to-stand functionality.



Back Recline



Leg Raise



Power Lumbar



Power Headrest



- **Contemporary Medium Size** available
- **Quattro Action:** Back recline, leg raise, lumbar support and headrest ideal for watching TV when reclined
- **Integrated Lift and Recline Functions** for easy sit-to-stand and maximum comfort
- **Soft Microsuede** upholstery finishes
- **Backlit Hand Control** with integrated USB charging point
- **Chaise Padding Design** fully supports legs, redistributing pressure



**2** YEAR WARRANTY



Latte Microsuede



Gunmetal Microsuede

| Code      | Colour   | Size   | Seat Height | Seat Depth | Seat Width | Backrest Height | Overall Height | Overall Width | SWL   |
|-----------|----------|--------|-------------|------------|------------|-----------------|----------------|---------------|-------|
| CHP227500 | Latte    | Medium | 480mm       | 530mm      | 530mm      | 630mm           | 1050mm         | 870mm         | 159kg |
| CHP227510 | Gunmetal | Medium | 480mm       | 530mm      | 530mm      | 630mm           | 1050mm         | 870mm         | 159kg |