




SIGNATURE 2 DUAL ACTION 2 MOTOR LIFT RECLINE CHAIR

The Signature 2 is a dual action lift recline chair offering independent movement of the backrest and legrest for optimal positioning. Available in three size options, the Signature 2 maximises comfort and support for a wide range of users.



*Custom Upholstery

- **Three sizes help match the chair's fit to the user's body shape**, improving comfort, posture and safety when standing
- **Dual action design creates infinite positioning** via independent backrest and legrest adjustment
- **Adjustable seat depth** ensures correct positioning and safer transfers for a wide range of users
- **Integrated leg extension** helps to comfortably support lower legs for a wide range of user heights
- **Swivel castor wheels** enable easy movement of the chair for cleaning or repositioning
- **All-day battery backup** ensures continued operation during power outages, supporting safety and independence
- **Backlit hand control and integrated USB charging** are simple and convenient

2 YEAR WARRANTY



Code	Colour	Size	Seat Height	Seat Depth	Seat Width	Backrest Height	Overall Height	Overall Width	SWL
CHP228200	Mink Warwick	Small	480mm	460 – 385mm	480mm	580mm	980mm	740mm	155kg
CHP228210	Lagoon Warwick								
CHP228220	Heritage								
CHP228230	Mink Warwick	Medium	480mm	490 – 515mm	510mm	620mm	1050mm	780mm	155kg
CHP228240	Lagoon Warwick								
CHP228250	Heritage								
CHP228260	Mink Warwick	Large	500mm	510 – 535mm	540mm	650mm	1110mm	820mm	155kg
CHP228270	Lagoon Warwick								
CHP228280	Heritage								

**Adjustable Neck Pillow CHP197781 Mink Warwick CHP197783 Lagoon Warwick CHP197782 Heritage