

MONTANA HD DUAL ACTION 2 MOTOR LIFT RECLINE CHAIR

This customisable bariatric lift recline chair offers dual-action positioning with 254kg or 450kg safe working load, a high-capacity legrest support, and configurable backrest options for specialised clinical applications and enhanced user seating experience.



Optional Backrests*



Immersion
Cushion



Waterfall
Cushion



Back Recline



Leg Raise



160kg Legrest
Weight Capacity**



Barley Vinyl



Custom Upholstery*



*Custom
Order Form

- **Dual Action:** Independent back recline and leg raise
- **True Bariatric Seating:** 254kg and 450kg options available
- **Custom Sizing:** Specify the seat width, depth, height and armrest height. This option is available via the custom order form*

- **Huge 160kg Legrest Weight Capacity:** loads of up to 160kg on HD254 model and 140kg on HD450 model
- **Configurable Backrests:** Optional alternative backrests are available on custom made chairs including laterally*



Code	Colour	Size	Seat Height	Seat Depth	Seat Width	Backrest Height	Overall Height	Overall Width	**Legrest SWL	SWL
CHP227405	Barley Vinyl	HD254	500mm	500mm	650mm	810mm	1270mm	900mm	160kg	254kg
CHP227410	Barley Vinyl	HD450	500mm	500mm	850mm	810mm	1270mm	1100mm	140kg	450kg

Note: Lifetime warranty frame, 3 year warranty upholstery and 2 year warranty mechanism and electronics