



# ODIN SCOOTER

## LARGE

**Comfortable Rides with Confidence, Wherever You Go**

The Odin Scooter delivers a smooth, stable, and secure ride, helping you enjoy your day with comfort and peace of mind. Designed for both daily errands and longer trips, it combines supportive seating, reliable suspension, and intuitive controls for a ride that feels effortless from start to finish.







### Reliable Performance

High-torque motor (700W) allows users to navigate steeper inclines.



### Plush, Supportive Seating

Features a deeply padded, adjustable captain's seat with suspension, and flexible tiller. Seat slides, swivels and reclines for perfect positioning.



### Smart Safety

Emergency stop button and automatic speed reduction while turning keeps you safe and secure.



### Comfortable Ride

Independent suspension absorbs bumps and reduces vibrations for comfort on uneven surfaces.



### Digital Dashboard

Clearly shows battery life, speed, and other key information, giving you confidence while you travel.



### Effortless Steering in Compact Spaces

A 1660mm turning circle allows smooth manoeuvring around shop aisles and community spaces.

Overall L x W x H	1420 x 680 x 1400mm
Seat Size D x W	440 x 530 - 580mm
Seat to Deck Height	470mm
Seat to Ground Height	670mm
Armrest Height	180 - 270mm
Min Turning Radius	1660mm
Ground Clearance	90mm
Max Climbing Angle	10°
Max Driving Speed	10km/h
Driving Range	40km
Tyres	330mm   Pneumatic
Battery	2 x 12V 75Ah (Lead Acid)
Product Weight (excl. battery)	92kg
Battery Weight	51.5kg
Power   Motor	700W   4 Pole
Controller	RHINO 2 120A
Brake System	Electro-Mechanical
Free Wheel Mode	Yes
Charger Size	24V 8A

\*Scooter Frame/Motors/Axles/Brakes 2 Years | Battery 6 Months | Electronics 1 Year



SCT671286  
Sky Blue



ASFY186.4.100925