



Oslo Seat Walker

Sleek lines, elegant curves, and eye-catching shades give the Oslo seat walker sophisticated poise. The aero-space grade aluminium frame is highly durable and ultra lightweight for easy manoeuvrability, folding, and lifting, while tool-free adjustments can suit individual user needs.

Min Weight
6.25kg



Ultra-lightweight Design: Constructed with aerospace-grade aluminium for exceptional strength, the 6.25kg minimum weight makes it easy to lift and transport

Simple to Adjust: Tool-free handle and backrest height adjustment accommodates users of varying stature, with push-button locking for secure customisation

Ergonomic Comfort: Contoured, non-slip hand grips to reduce wrist strain and adjustable backrest for tailored support

All-terrain Mobility: Puncture-resistant PU wheels with swivel front castors and fixed rear wheels improve manoeuvrability

Portable and Compact: One-click folding mechanism collapses the rollator into a slim profile (20cm width) for effortless car boot storage or travel

User-centric Accessories: Secure zippered storage bag with 5kg capacity

	
Purple WAF705320PU	Green WAF705320GR

2 YEAR WARRANTY
Aluminium Frame

1 YEAR WARRANTY
Components



Overall Height mm	Overall Length mm	Overall Width mm	Seat Width mm	Seat Height mm	Seat Depth mm	Folded Width mm	Front Wheel mm	Rear Wheel mm	Min Weight kg	Total Weight kg
830 - 940	660	540	450	585	205	200	230	195	6.25	6.9

1300 133 120 | aspirecare.com.au