



PostureFit Lift Recline Chair

Configurable Pressure Care Seating Solution for unique and changing needs

The Aspire PostureFit is a configurable powerlift and recline chair with multiple options for a therapist or carer to ensure the best fit for a client's postural and pressure care needs. The PostureFit is available in two colours and three sizes to suit a range of users.



Zero Gravity



Back Recline



Leg Raise



Power Headrest



Vertical-lift Mechanism

CARRFLEX Available in Black and Stone medical grade multistretch upholstery

Tri-action Chair: Features powered headrest in addition to the independent back and leg motion

Zero-gravity Positioning: Relieves lower back pressure for pain management when fully reclined

Integrated Footrest Extension: Fits most user heights without altering the chair's legrest or height

Tool-free Seat Depth and Width Adjustment: Easily adjust the backrest for two seat depths and interchange armrests for two seat widths, both with a 50mm range

Medical Grade Foam: High density foam allows for prolonged periods of sitting

Backrest Cushions: Have adjustable filling and support wedges are also available to meet unique postural requirements

Effortless Mobility and Safety: Large castors engage via a rear foot pedal, ensuring stability and preventing tipping

Transport Mode: Electrical functions deactivate for safety

Battery Backup: To maximise safety and independence



PostureFit Brochure
view optional accessories

2 YEAR WARRANTY
Frame and Electronics

1 YEAR WARRANTY
Handset and Cables



Code	Description	Overall Depth	Overall Width	Overall Height	Seat Depth	Seat Width	Seat Height*	Backrest Height
 CHP198240 CHP198241	Petite	900mm	760 - 810mm	1080mm	475 - 525mm	485 - 535mm	450mm	700mm
 CHP198250 CHP198251	Small/Medium	900mm	760 - 810mm	1140mm	475 - 525mm	485 - 535mm	490mm	740mm
 CHP198260 CHP198261	Medium/Large Tall	950mm	820 - 870mm	1180mm	525 - 575mm	535 - 585mm	490mm	790mm

1300 133 120 | aspirecare.com.au

*Seat-to-floor Height