



# Raphael X Quattro Action

## 4 Motor Lift Recline Chair

Designed for enhanced comfort and customisable support. The Raphael X delivers exceptional comfort with features to suit a diverse range of users. The adjustable backrest includes zipper access for personalised postural adjustment and lateral support, providing improved upright sitting. It's plush design compliments any modern home environment.



**Lateral Backrest**  
Standard



**Classic Backrest**  
Optional\*

### Colour range



Golden Grain  
Microsuede



Marina Grey  
Microsuede



Nordic  
Stone Vinyl



Custom  
Upholstery\*

### Functions



Back  
Recline



Leg  
Raise



Power  
Lumbar



Power  
Headrest



\*Custom Order Form

The Raphael X range is available in five sizes—Petite, Small, Medium, Large, and Maxi. Each size is designed to provide the right balance of comfort, postural support, and ease of use, ensuring a well-fitted seating experience for a wide range of users.



**Quattro Action Mechanism**

Independently adjusts the backrest, leg rest, lumbar support, and headrest to provide tailored comfort for various activities, including reclining, relaxing, and seated support

**Adjustable Backrest for Personalised Support:** Backrest cushions can be repositioned to accommodate individual postural needs, promoting better alignment and comfort

**Built-in Lateral Supports:** Helps maintain an upright sitting posture, improving comfort and stability

**Integrated Footrest Extension:** Helps with fatigue and improves lower limb comfort<sup>1</sup>

**Adjustable Seat Depth<sup>2</sup>:** Offers 25mm of depth adjustment, allowing taller users to achieve a better-seated position

**Protective Covers for High-Wear Areas:** Removable covers for the headrest and armrests to keep your chair looking new



**Backlit Hand Control with USB Charging**

Easy-to-use, illuminated controls improve visibility in low light, while the integrated USB port allows convenient device charging without needing to leave the chair

**Refined, Contemporary Design:** Larger, plush armrests provide extra comfort while maintaining a streamlined look that suits home environment

**All-day Battery Backup:** A lithium-ion battery ensures continued operation during power outages, supporting safety and independence

**Premium Upholstery Options:** Available in modern finishes, including microsuede and leather-look vinyl, with local customisation available by request

**Classic Backrest Option:** Available as an additional accessory for users who prefer a more traditional support profile



2

YEAR

WARRANTY



SWL  
130kg



MAXI SWL  
180kg

Code	Colour	Size	Seat Height mm	Seat Depth <sup>2</sup> mm	Seat Width mm	Seat to Top of Back mm	Overall Height mm	Overall Depth mm	Overall Width mm	Distance from Wall mm	SWL kg
CHP227640	Golden Grain	Petite*	450	450-475	470	630	1060	840	820	420	130
CHP227641	Marina Grey										
CHP227642	Nordic Stone										
CHP227650	Golden Grain	Small	450	480-505	500	630	1060	870	850	400	130
CHP227651	Marina Grey										
CHP227652	Nordic Stone										
CHP227660	Golden Grain	Medium	470	500-525	525	680	1100	900	875	400	130
CHP227661	Marina Grey										
CHP227662	Nordic Stone										
CHP227670	Golden Grain	Large	500	530-555	550	720	1100	910	900	420	130
CHP227671	Marina Grey										
CHP227672	Nordic Stone										
CHP227680	Golden Grain	Maxi	510	540-565	650	675	1120	910	1020	420	180
CHP227681	Marina Grey										
CHP227682	Nordic Stone										

<sup>1</sup>The Petite model is designed without a leg extension, catering to the needs of shorter users.  
<sup>2</sup>Only available on stock produced from February 2025 onwards