



# Raphael X5

## Five Action Lift Recline Chair

The Raphael X5 lift recline chair features individual control of all components allowing you to select the perfect comfort position to suit your needs. With customisable material options available, the Raphael X5 is guaranteed to compliment your personal style.



*Shown in Marina Grey Microsuede*



Zero Gravity  
Tilt-in-space



Back  
Recline



Leg  
Raise



Power  
Lumbar



Power  
Headrest



**Zero Gravity and Tilt-in-space:** Create a feeling of weightlessness by elevating your legs above your upper body, with only a limited amount of space needed behind the lift recline chair



**Legrest Extension:** Provides improved support under the lower legs, especially beneficial for taller users<sup>1</sup>



**Electric Adjustable Headrest:** Ideal for watching TV, reading, or relaxing while reclined, ensuring your neck is always well supported



**Variable Seat Depth:** A 25mm range of seat depth ensures users achieve an optimal seating position



Simplified Handset

**Precise and Simple Handset Options:** Achieve the perfect position with the precise handset. An optional simplified handset is available for basic functions. Handset loops allow for easy access on either side of the chair





**Custom Upholstery**

Choose from a range of leather, microsuede, and fabric options for a custom finish and personalised style



View our full range



Golden Grain Microsuede



Marina Grey Microsuede



Nordic Stone Vinyl



Ebony Leather

**SIZING CHART**

Size	Seat Height	Seat Depth	Seat Width	Armrest Height	Backrest Height	Overall Height	Overall Depth	Overall Width	Distance from Wall
Petite <sup>1</sup>	450mm	450-475mm	470mm	160mm	630mm	1060mm	840mm	820mm	420mm
Small	450mm	480-505mm	500mm	175mm	630mm	1060mm	870mm	850mm	400mm
Medium	470mm	500-525mm	525mm	175mm	680mm	1100mm	900mm	875mm	400mm
Large	500mm	530-555mm	550mm	175mm	720mm	1160mm	910mm	900mm	420mm
Maxi	510mm	540-565mm	650mm	175mm	675mm	1120mm	900mm	1020mm	420mm

**CODE CHART**

Colour	Petite	Small	Medium	Large	Maxi
Golden Grain	CHP228400	CHP228500	CHP228600	CHP228700	CHP228800
Marina Grey	CHP228410	CHP228510	CHP228610	CHP228710	CHP228810
Nordic Stone	CHP228420	CHP228520	CHP228620	CHP228720	CHP228820
Ebony	—	CHP228530	CHP228630	CHP228730	—

**ACCESSORIES**

Optional Simplified Handset CHP228900

<sup>1</sup>The Petite model is designed without a leg extension, catering to the needs of shorter users