



Vogue Forearm Seat Walker

Lightweight Aluminium Body Improves manoeuvrability and portability

For those who need a little more support, the Forearm Walker is for you. Soft and contoured forearm pads ensure maximum comfort and a smooth ride.



Smart Posture System
With adjustable forearm supports across height, depth and angle



Simple Folding Mechanism
Folds to 225mm



Removable Storage
Zip bag secures up to 5kg of essentials

Luxury Look and Feel Danish design, profiled tubing, moulded fittings, ergonomic handles, and refined colours

Padded Forearm Support For enhanced comfort

Lightweight and Ergonomic With small turning circle for excellent manoeuvrability

Advanced Locking Handbrakes For increased user safety when walking and sitting

AutoLocate Smooth Height Adjustment That automatically locates to pre-set heights, making adjustment simple and accurate

Unbreakable TPE Super-soft Wheels Offers high durability and a smooth ride for increased comfort

Bonus Walking cane holder

1 YEAR WARRANTY

SWL 150kg

SWL 23.6st

White / Black
WAF705400

Overall Length mm / "	Overall Width mm / "	Overall Height mm / "	Seat Width mm / "	Seat Height mm / "	Folded Width mm / "	Wheel Size Front & Rear mm / "	Minimum Functional Weight kg / lbs	Typical Functional Weight* kg / lbs
655	620	900 - 1170	441	620	225	200	9.2	10
26	24	40 - 50	17	24.5	9	8	20	22

AS/NZ Approved | Vogue is registered in the ARTG, AUST R ARTG 271255 | *Includes accessories, bag, and backstrap.

Australia
1300 133 120 | aspirecare.com.au

United Kingdom
01622 541235 | aidacare.co.uk