



Vogue Lightweight 2 Seat Walker

Lightweight Aluminium Body Improves manoeuvrability and portability

This modern European designed seat walker has a clean and crisp frame design and is available in a selection of eye-catching colours and prints to really stand-out while keeping you on the move.


Less Than
6kg



Handbrakes Lock
For user safety and stability



Compact Turning Circle
For improved mobility in tight spaces



Leopard
WAF705310LE



Flower
WAF705310FL



Zebra
WAF705310ZB



Red
WAF705310RE



Champagne
WAF705310CH



Ergonomic Handles With personalised height adjustment to promote comfortable and correct posture

Deep Seat and Wide Backrest For improved comfort and support
Bonus Walking cane holder



Adjustable Handles
Allows for personalised positioning, helping to maintain proper posture



Simple Folding Mechanism
Locks for easy transportation and storage



Removable Storage
Zip bag secures up to 5kg of essentials



Accessories
With a cane holder included as standard, choose from our range of optional accessories

Overall Length mm	Overall Width mm	Overall Height mm	Seat Depth mm	Seat Width mm	Seat Height mm	Folded Width mm	Wheel Size mm	Minimum Functional Weight kg	Typical Functional Weight* kg
655	620	780 - 920	240	465	540	220	200	5.4	6.45

AS/NZ Approved | Vogue is registered in the ARTG, AUST R ARTG 271255 | *Includes accessories, bag, and backstrap.