

ComfiMotion Collection

Luxurious adjustable bed collection for an enhanced sleeping experience



 aspire ComfiMotion®

ComfiMotion

collection

We believe that creating a supportive and luxurious sleep environment can significantly improve health and wellbeing.

Our ComfiMotion Collection is designed to blend seamlessly into the home environment, whilst improving comfort and independence.


Create Your ComfiMotion



Bed Function



Headboard Style



Mattress Support



Accessories

© Aidacare Pty Ltd 2025.
No part of this brochure may be copied, distributed, reproduced, scanned or stored on any electronic database, whether in whole or in part in any form or by any means.


Due care has been taken to ensure all the information contained in this brochure is correct at the time of printing. All measurements, pictures, colours and weight capacities are to be used as a guide only. We reserve the right to modify the design or appearance of any product displayed in this brochure without prior notice.

ComfiMotion is a registered trademark of Aidacare Pty Ltd.

Activ Care

Ultra Low Adjustable Bed

Function & Style Support



Supporting the modern, independent lifestyle.

The Activ Care bed combines the best of in-home care, functionality, positioning and low height safety with luxurious and contemporary style for enhanced comfort, care and wellbeing.

Features

- **Ultra Low Height:** Minimum platform height of 195mm reduces falls risk, improving user safety.
- **Customisable Safe Entry/Exit Height:** Programmable single button press ensures perfect bed height to maximise independence.
- **Underbed Clearance:** Provides clear access for patient lifters and overbed tables.
- **Battery Backup:** Provides confidence and maintains independence during power outages.
- **Night Assist:** Luminous underbed halo for safe entry and exit.
- **Illuminated Handset:** Large button interface is intuitive and highly visible at all times.
- **Massage:** Relaxing vibration modes for personalised neck, shoulder and back relief.

Functions

- Hi-Lo Functionality:** Assists with manual handling and helps to mitigate fall risks.
- Head Tilt:** Relaxes neck muscles, providing optimal support whilst reading or watching TV.
- Back Raise:** Improves postural support and can ease symptoms of reflux, respiratory conditions, snoring and sleep apnoea.
- Auto-Regression in the Backrest:** Stops the user sliding down the bed during movement while improving comfort and maximising pressure redistribution.
- Clinical Positioning:** One touch raises back and legs to alleviate cardiac exertion, improving resting comfort.

- Leg Raise:** Promotes circulation, relieving tired, swollen, fatigued or painful limbs.
- Memory Function:** Allows the user to set and return to their preferred position at any time with one click of a button.
- Zero Gravity:** Helps to relax the circulatory system and muscles of the lower back.
- Seated Position:** Assists with bed entry and exit and upright dining.
- Trendelenburg (and Reverse):** Can assist care givers when repositioning a user*.
- Lockable Castors:** Eight heavy duty dual castors.

Upholstery



Portobello Grey

Sizes (Length x Width x Height Range)

Long Single 2030 x 920 x 195 – 700mm BEB046740	Long Double 2030 x 1370 x 195 – 700mm BEB046760
King Single 2030 x 1070 x 195 – 700mm BEB046750	Queen 2030 x 1520 x 195 – 700mm BEB046770



▲ Companion Sleeping
See Configurations section for more companion options.

*Carers should seek advice from a trained Allied Health professional for correct manual handling techniques.

10 YEAR LIMITED WARRANTY

SWL
250kg
LS & KS


SWL
360kg
LD & Q

Visit comfimotion.com.au for Limited Warranty T&Cs

Care

Adjustable Bed

Function & Style Support



For comfort you deserve.

The Care adjustable bed blends seamlessly into the home without compromising on the therapeutic needs of the user or safety of partners or carers. Select from multiple position options designed to maximise independent mobility and improve sleep quality.

Features

- **Night Assist:** Luminous underbed halo for safe entry and exit throughout the night.
- **Illuminated Handset:** Large button interface is intuitive and highly visible at all times.
- **Enclosed Frame Design:** Ensures optimal positioning when getting into bed and safely covers the moving base platform.
- **Safe-Lowering Technology:** Ensures that back, leg and Hi-Lo functions only lower under the force of gravity, reducing potential for accidents or injuries.
- **Box-on-Box Design:** Facilitates easier entry and exit and allows true flush positioning of beds in companion or dual adjustable bed configurations.

Functions

Hi-Lo Functionality: Assists with manual handling and helps to mitigate fall risks.

Back Raise: Improves postural support and can ease symptoms of reflux, respiratory conditions, snoring and sleep apnoea.

Clinical Positioning: Raises back and legs to alleviate cardiac exertion, improving resting comfort.

Lockable Castors: Four heavy duty castors.


Leg Raise: Promotes circulation, relieving tired, swollen, fatigued or painful limbs.

Zero Gravity: Helps to relax the circulatory system and muscles of the lower back.

Seated Position: Assists with bed entry and exit and upright dining.

Trendelenburg (and Reverse): Can assist care givers when repositioning a user*.


Upholstery



Portobello Grey

*Carers should seek advice from a trained Allied Health professional for correct manual handling techniques.

Variable Bed Height Options



Low Height
280 – 580mm
with Low Profile
Bed Legs

Standard
350 – 650mm
with Standard
Castors

Accessibility
400 – 700mm
with Large
Castors (optional)



▲ Companion Sleeping

See Configurations section for more companion options.

Sizes (Length x Width x Height Range)

Long Single – Fixed Companion Bed Available
2030 x 920 x 280 – 580mm¹ BEB046720

King Single
2030 x 1070 x 280 – 580mm¹ BEB046730

1. With Low Profile Bed Legs.



Visit comfimotion.com.au for Limited Warranty T&Cs

Relax

Fixed Height Adjustable Bed

Function & Style Support



Supporting the modern, independent lifestyle.

The Relax Bed has a fixed base and combines luxurious finishes with multiple positioning and therapeutic functions to support users and carers alike.

Features

- **Night Assist:** Luminous underbed halo for safe entry and exit throughout the night.
- **Massage:** Relaxing vibration modes for personalised neck, shoulder and back relief.
- **Press & Hold Safety Handset:** Intuitive, large button interface enhances usability.
- **Box-on-Box Design:** Facilitates easier entry and exit, allowing true flush positioning of beds in companion or dual adjustable bed configurations.
- **German Engineered Motors:** For smooth, quiet adjustment.
- **Personalised Positioning:** Hand control allows users to find their perfect position to maximise comfort and relaxation.
- **USB Charge:** Connect and charge your devices conveniently.

Functions

Head Tilt: Relaxes neck muscles, providing optimal support whilst reading or watching TV.

Back Raise: Improves postural support and can ease symptoms of reflux, respiratory conditions, snoring and sleep apnoea.

Leg Raise: Promotes circulation, relieving tired, swollen and/or painful limbs.

Zero Gravity: Helps to relax the circulatory system and muscles of the lower back.

Memory Function: Allows the user to set and return to their preferred positions at any time with one click of a button.

Upholstery



Portobello Grey

Sizes (Length x Width x Height Range)

Long Single
2030 x 920 x 280 – 700mm **BEB046900**


King Single
2030 x 1070 x 280 – 700mm **BEB046910**

Long Double
2030 x 1370 x 280 – 700mm **BEB046920**

Long Queen
2030 x 1520 x 280 – 700mm **BEB046930**



▲ Companion Sleeping
See Configurations section for more companion options.



Visit comfimation.com.au for Limited Warranty T&Cs

Mattresses

Engineered for Motion

Our ComfiMotion Pocket Spring and Memory Foam mattresses are designed with the flexibility to follow the contouring profile of an adjustable bed.

Pocket Spring

The ComfiMotion Pocket Spring mattress offers zoned support springs that respond and contour to the individual body shape to improve comfort, maximising support and sleep quality.




Transitional Layered Pocket Springs arranged in weight activated support zones.



Integrated Side-Walls create roll resistance at the bed edge and provide extra rigidity to ease bed entry and exit.



Luxuriously Soft Pillowtop Layer provides a plush immersive feeling complimented by the responsive pocket springs.



▲ Pocket Spring Mattress


▲ Classic Memory Foam Mattress



Contoured Immersion: Specially layered memory foam provides contoured immersion for support and pressure relief.



Superior Breathability: Up to 8x more breathable than conventional memory foams.



Anti-microbial: Environmentally safe treatment is added to the foam which reduces odours and kills 99.9% of bacteria.



Pressure Care Incontinence Cover
Breathable, 4-way stretch waterproof cover to reduce shear, friction and maintain skin micro-climate.



Integrated Side-Walls create roll resistance at the bed edge and provide extra rigidity to ease bed entry and exit.

Pocket Spring


(Length x Width x Height)

Pocket Spring Mattress – Long Single
2030 x 920 x 240mm [BEM469800](#)

Pocket Spring Mattress – King Single
2030 x 1070 x 240mm [BEM469810](#)

Pocket Spring Mattress – Long Double
2030 x 1370 x 240mm [BEM469820](#)

Pocket Spring Mattress – Queen
2030 x 1520 x 240mm [BEM469830](#)



Visit comfimotion.com.au for Limited Warranty T&Cs

Memory Foam - Classic

(Length x Width x Height)

Long Single Classic
2030 x 920 x 200mm [BEM469910](#)

King Single Classic
2030 x 1070 x 200mm [BEM469920](#)

Long Double Classic
2030 x 1370 x 250mm [BEM469922](#)

Queen Classic
2030 x 1520 x 250mm [BEM469924](#)

Memory Foam - Pressure Care

(Length x Width x Height)

Long Single Pressure Care
2030 x 920 x 200mm [BEM469930](#)

King Single Pressure Care
2030 x 1070 x 200mm [BEM469940](#)

Long Double Pressure Care
2030 x 1370 x 250mm [BEM469942](#)

Queen Pressure Care
2030 x 1520 x 250mm [BEM469944](#)

PRESSURE CARE



MEDIUM



Visit comfimotion.com.au for Limited Warranty T&Cs

Clinical Considerations

Accessibility and Inclusivity: Unexposed bed surround with padded upholstery removes gaps between beds in co-sleeping arrangements and promotes safe entry and exit.

User Independence: Height adjustability improves user mobility and assists carers with manual handling.

Flexible Positioning: Caters to various user requirements such as respiratory conditions, snoring, sleep apnoea, pressure redistribution, joint or muscle pain, lymphoedema and swelling.

Vision and Hand Function: Handset with a clear and intuitive button display, assists users with reduced vision or hand function.

Falls Prevention: Integrated bed lighting and a minimum platform height of 195mm (Activ Care) and 280mm (Care and Relax) helps to reduce the likelihood and severity of falls in the home.

User Safety: Battery backup provides peace of mind in the event of a power outage (Standard on Activ Care, optional accessory on Care).





Clinical Assistance Guide

To view our detailed clinical guide, please visit aidacare.com.au or scan the QR code


Headboards



Allure

Stylish freestanding headboard
Subtle square padded design

Long Single **BEA830120** Queen **BEA830134**
King Single **BEA830130** King **BEA830136**
Long Double **BEA830132**



Chesterfield

Stylish freestanding headboard
Chesterfield padded design

Long Single **BEA830100** Queen **BEA830114**
King Single **BEA830110** King **BEA830116**
Long Double **BEA830112**

Pillows



Breeze Cooling Gel Memory Foam

- Adaptive and supportive sleep surface
- Helps regulate sleep surface temperature
- Accommodates side and back sleepers

BEP050100



Plush Memory Foam and Fibrefill

- Familiar sleep surface
- Improved head and neck support
- Accommodates side and back sleepers

BEP050110

Accessories


Function & Style Support



Crook Handle Bed Stick

Care/Relax Bed **BEA830200**
Activ Care Bed **BEA830211**


Provides multiple grip points to assist repositioning in bed, and safe entry and exit.



Extendable Side Rail

Activ Care Bed **BEA830205**


Provides a telescoping bed rail support for users at risk of falls.



Bed Rail

Care/Relax Bed **BEA830210**
Activ Care Bed **BEA830212**

Improves bed edge awareness and offers multiple grip points for repositioning and bed mobility.



Split King Upgrade Kit

Activ Care Bed LS & KS only **BEA830230**

Ensures simultaneous Hi-Lo platform movement whilst still allowing individual positioning requirements when beds are in partner configuration.



Simplified Handset

Activ Care Bed **BEA830214**


For clients with reduced cognitive ability.



Partner Handset

Activ Care Bed **BEA830213**

Allows double and queen beds to be operated from either side of the bed.



Standard




Large

Large Accessibility Castors

Care Bed only (4 x required) **MWT451080**

Larger castors raise the bed frame to 155mm, facilitating overbed table and patient lifter access.



Battery Backup

Care & Relax beds (shown) **BEA830220**
Activ Care beds **BEB046590**

Provides peace of mind in the event of a power outage.

Collection Summary


Beds

Care Adjustable Bed – (LS)	BEB046720	Activ Care Adjustable Bed – (LS)	BEB046740
Care Adjustable Bed – (KS)	BEB046730	Activ Care Adjustable Bed – (KS)	BEB046750
Relax Adjustable Bed – (LS)	BEB046900	Activ Care Adjustable Bed – (LD)	BEB046760
Relax Adjustable Bed – (KS)	BEB046910	Activ Care Adjustable Bed – (Q)	BEB046770
Relax Adjustable Bed – (LD)	BEB046920	Fixed Companion Bed – (LS)	BEB046800
Relax Adjustable Bed – (Q)	BEB046930		



Mattresses

Pocket Spring Mattress – (LS)	BEM469800	Memory Foam Mattress – Pressure Care (LS)	BEM469930
Pocket Spring Mattress – (KS)	BEM469810	Memory Foam Mattress – Pressure Care (KS)	BEM469940
Pocket Spring Mattress – (LD)	BEM469820	Memory Foam Mattress – Pressure Care (LD)	BEM469942
Pocket Spring Mattress – (Q)	BEM469830	Memory Foam Mattress – Pressure Care (Q)	BEM469944
Memory Foam Mattress – Classic (LS)	BEM469910		
Memory Foam Mattress – Classic (KS)	BEM469920		
Memory Foam Mattress – Classic (LD)	BEM469922		
Memory Foam Mattress – Classic (Q)	BEM469924		



Headboards

Allure Headboard – (LS)	BEA830120	Chesterfield Headboard – (LS)	BEA830100
Allure Headboard – (KS)	BEA830130	Chesterfield Headboard – (KS)	BEA830110
Allure Headboard – (LD)	BEA830132	Chesterfield Headboard – (LD)	BEA830112
Allure Headboard – (Q)	BEA830134	Chesterfield Headboard – (Q)	BEA830114
Allure Headboard – (K)	BEA830136	Chesterfield Headboard – (K)	BEA830116



Pillows

Breeze Pillow	BEP050100	Plush Pillow	BEP050110
---------------	---------------------------	--------------	---------------------------



Accessories

Crook Handle Bed Stick – Care/Relax Bed	BEA830200	Large Accessibility Castors – Care Bed	MWT451080
Crook Handle Bed Stick – Activ Care Bed	BEA830211	Battery Backup – Care/Relax Bed	BEA830220
Bed Rail – Care/Relax Bed	BEA830210	Split King Upgrade Kit – Activ Care Bed	BEA830230
Bed Rail – Activ Care Bed	BEA830212	Simplified Handset – Activ Care Bed	BEA830214
Extendable Bed Rail	BEA830205	Partner Handset – Activ Care Bed	BEA830213



Configurations for Companions


Whichever ComfiMotion package you select, rest assured there is a companion solution available to suit those wishing to sleep alongside any ComfiMotion Adjustable bed.

What to Order


KING
2030 x 1840



KING PLUS
2030 x 1990



SUPERKING
2030 x 2140



Combine any two singles below:

Care Long Single	BEB046720
Active Care Long Single	BEB046740
Fixed Long Single	BEB046800
Relax Long Single	BEB046900


Combine any single with king single below:

Relax Long Single	BEB046900
Care Long Single	BEB046720
Active Care Long Single	BEB046740
Fixed Long Single	BEB046800
Relax King Single	BEB046910
Care King Single	BEB046730
Active Care King Single	BEB046750

Combine any two king singles below:


Care King Single	BEB046730
Active Care King Single	BEB046750
Relax King Single	BEB046910

Long Double and Queen Beds can not be paired with other beds in the range.



Dual Adjustable Beds

Designed to sit flush when positioned together to enable close monitoring and enhanced care. This configuration offers the flexibility of independent adjustability, without compromising on the sense of inclusivity and intimacy between partners.



Split King Upgrade Kit

Ensures simultaneous Hi–Lo platform movement whilst still allowing individual positioning requirements when beds are in partner configuration.



Fixed Companion Bed

Matching single bed without any adjustability designed to sit flush when positioned together. It facilitates a side–by–side sleeping arrangement at a fixed height of 370mm, without compromising on the users care and assistance requirements.

Split King Upgrade Kit
[BEA830230](#)

Compatible with active care bed

Fixed Companion Bed – (LS)
2030 x 920 x 370mm [BEB046800](#)

LS – Long Single

KS – King Single

LD – Long Double

Q – Queen

K – King

AIDACARE PTY LTD

Building 3A, 1 Moorebank Avenue,
Moorebank NSW 2170
Australia

1300 133 120 | aidacare.com.au

AIDACARE LIMITED

Arcadia House, Maritime Walk,
Ocean Village, Southampton SO14 3TL
United Kingdom

01622 541235 | aidacare.co.uk

MADE IN PRC

ASBR28.101 160725

