

aspire™ LIFT RECLINE CHAIRS



	GOLD RANGE															
	Gabriel X			Raphael X				Rembrandt		Matisse		Da Vinci		Renoir		Monet
Functions	Five Action			Quattro Action				Five Action		Quattro Action		Quattro Action		Quattro Action		Dual Action
Motors	5 Motor			4 Motor				5 Motor		4 Motor		4 Motor		4 Motor		3 Motor
Codes (Bold numbers refer to upholstery)	19 -CHP227590	20 -CHP227595	1 -CHP227640 2 -CHP227641 3 -CHP227642	1 -CHP227650 2 -CHP227651 3 -CHP227652	1 -CHP227660 2 -CHP227661 3 -CHP227662 21 -CHP227663	1 -CHP227670 2 -CHP227671 3 -CHP227672	1 -CHP227680 2 -CHP227681 3 -CHP227682	6 -CHP227550	7 -CHP227555	7 -CHP227560	4 -CHP227500 5 -CHP227510	8 -CHP227575	8 -CHP227580	8 -CHP227585	9 -CHP227540	
Size	Small	Large	Petite	Small	Medium	Large	Maxi	Large	Petite	Large	Medium	Small	Medium	Large	Large	
Seat Height	475mm	500mm	450mm	450mm	470mm	500mm	510mm	495mm	455mm	500mm	480mm	450mm	470mm	470mm	495mm	
Seat Depth	450mm	520mm	450 - 475mm	480 - 505mm	500 - 525mm	530 - 555mm	540 - 565mm	535mm	455mm	546mm	500mm	510mm	530mm	530mm	535mm	
Seat Width	490mm	530mm	460mm	480mm	500mm	530mm	540mm	560mm	510mm	571mm	530mm	510mm	550mm	590mm	560mm	
Backrest Height	710mm	790mm	630mm	630mm	680mm	720mm	675mm	600mm	560mm	730mm	630mm	660mm	760mm	760mm	570mm	
Overall Height	1100mm	1170mm	1100mm	1170mm	1100mm	1100mm	1120mm	1095mm	1015mm	1067mm	1050mm	1010mm	1060mm	1060mm	1065mm	
Overall Width	850mm	900mm	820mm	850mm	875mm	900mm	1020mm	1015mm	915mm	965mm	870mm	750mm	800mm	830mm	890mm	
SWL	150kg	155kg	130kg	130kg	130kg	130kg	180kg	159kg	136kg	159kg	159kg	130kg	130kg	130kg	159kg	
Warranty	2 Years	2 Years	2 Years	2 Years	2 Years	2 Years	2 Years	2 - 5 Years*	2 Years	2 Years	2 Years	2 Years	2 Years	2 Years	2 - 5 Years*	



	SILVER RANGE						BRONZE RANGE						
Functions	Signature 2			Signature 2			Montana X Maxi		Montana HD		Idaho	Oregon	Maui
	Single Action			Dual Action			Single Action		Dual Action		Dual Action	Single Action	Dual Action
Motors	Space Saver 1 Motor			2 Motor			1 Motor		2 Motor		2 Motor	1 Motor	2 Motor
Codes (Bold numbers refer to upholstery)	10-CHP228100	10-CHP228130	10-CHP228160	10-CHP228200	10-CHP228230	10-CHP228260	1-CHP227401	18-CHP227405	18-CHP227410	16-CHP227160 17-CHP227165	16-CHP227150 17-CHP227155	14-CHP227110	
	11-CHP228110	11-CHP228140	11-CHP228170	11-CHP228210	11-CHP228240	11-CHP228270							
	12-CHP228120	12-CHP228150	12-CHP228180	12-CHP228220	12-CHP228280	12-CHP228280							
Size	Small	Medium	Large	Small	Medium	Large	Maxi	HD254	HD450	Medium	Medium	Medium	
Seat Height	480mm	480mm	500mm	480mm	480mm	500mm	510mm	500mm	500mm	480mm	480mm	500mm	
Seat Depth	460 - 485mm	490 - 515mm	510 - 535mm	460 - 485mm	490 - 515mm	510 - 535mm	530mm	500mm	500mm	480mm	480mm	500mm	
Seat Width	480mm	510mm	540mm	480mm	510mm	540mm	640mm	650mm	850mm	480mm	480mm	540mm	
Backrest Height	580mm	620mm	650mm	580mm	620mm	650mm	670mm	810mm	1100mm	620mm	620mm	630mm	
Overall Height	980mm	1050mm	1110mm	980mm	1050mm	1110mm	1060mm	1270mm	1270mm	1100mm	1100mm	1100mm	
Overall Width	740mm	780mm	780mm	740mm	780mm	780mm	940mm	900mm	1100mm	800mm	800mm	780mm	
SWL	145kg	145kg	145kg	145kg	155kg	155kg	220kg	254kg	450kg	150kg	150kg	130kg	
Warranty	2 Years	2 Years	2 Years	2 Years	2 Years	2 Years	2 Years	2 - 3 Years		2 Years	2 Years	2 Years	

Tested to International Standards | Enquire for DVA Contracted Lift Recline Chairs | *2 year components and 5 year frame warranty | 2 year limited warranty, for details refer to respective lift recline chair user manuals.

Upholstery Range

1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21
Golden Grain Microsuede	Marina Grey Microsuede	Nordic Stone Vinyl	Latte Microsuede	Gunmetal Microsuede	Canyon Walnut Microsuede	Santa Barbara Fog Microsuede	Dusk Microsuede	Canyon Walnut Microsuede	Mink Suede	Lagoon Suede	Heritage Vinyl	Fawn Soft Suede	Mocha Chenille Fabric	Mushroom Chenille Fabric	Charcoal Fabric	Latte Vinyl	Barley Vinyl	Dark Grey	Ebony	Stone Leather