



Vogue Adventure Seat Walker

Equipped with Extra-large All-terrain Wheels Built for travel beyond the everyday

Rugged and durable, the Adventure Seat Walker is the perfect all terrain walker. Super-sized front wheels can manage any surface, perfect for the adventurous user.



Large Multi-terrain Wheels
For durability and stability on and off road



Simple Folding Mechanism
Locks for easy transportation and storage



Removable Storage
Zip bag secures up to 5kg of essentials

Ergonomic Design Offers personalised height adjustment which promotes correct posture

Heavy-duty Frame and Wider Seat Provide enhanced support for users up to 200kg

Handbrakes Lock For user safety and stability

Compact Turning Circle For improved mobility in tight spaces

Bonus Walking cane holder

1 YEAR WARRANTY



Black
WAF705450

Overall Length mm	Overall Width mm	Overall Height mm	Seat Width mm	Seat Height mm	Folded Width mm	Wheel Size Front & Rear mm	Minimum Functional Weight kg	Typical Functional Weight* kg
770	695	900 - 1035	540	530	245	279 & 225	8.8	10

AS/NZ Approved | Vogue is registered in the ARTG, AUST R ARTG 271255 | *Includes accessories, bag, and backstrap.

1300 133 120 | aspirecare.com.au