



Vogue Seat Walker Range

Premium Mobility Solutions



VOGUE CARBON FIBRE SEAT WALKER

Carbon Fibre Frame Extremely lightweight to maximise portability

One of the lightest seat walkers on the market. Our Vogue Carbon Fibre Seat Walker epitomises function and style with a super-lightweight carbon fibre frame for a superior walking experience.

Less Than
5kg



Compact Turning Circle
For improved mobility in tight spaces



Deep Seat and Wide Backrest
For improved comfort and support



Handbrakes Lock
For user safety and stability

Ergonomic Design Offers personalised height adjustment which promotes correct posture

Removable Storage Zip bag Secures up to 5kg of essentials

Simple Folding Mechanism Locks for easy transportation and storage

Bonus Walking cane holder

SWL
150kg

SWL
23.6st

1 YEAR WARRANTY


Sapphire
Short


Black
Medium


Emerald Green
Tall

	Overall Length	Overall Width	Overall Height	Seat Width	Seat Height	Folded Width	Wheel Size	Minimum Functional Weight	Typical Functional Weight**
	mm / "	mm / "	mm / "	mm / "	mm / "	mm / "	mm / "	kg / lbs	kg / lbs
Sapphire - Short WAF705340	655	600	715 - 830	460	500	230	203	4.7	5.9
	26	24	28 - 32	18	20	9	8	10.4	13
Black - Medium WAF705350	655	615	725 - 940	460	550	230	203	4.9	6.2
	26	24	28 - 37	18	22	9	8	11	13.2
Emerald Green - Tall WAF705370	655	600	780 - 1030	460	620	230	203	5.2	6.3
	26	24	31 - 41	18	24	9	8	11.5	14

VOGUE NINJA CARBON FIBRE SEAT WALKER

Carbon Fibre Frame Extremely lightweight to maximise portability

The Vogue Ninja is a lightweight seat walker designed for everyday use. Made from carbon fibre, it's easy to lift in and out of a car and has large front wheels for greater mobility outdoors.

Less Than
5.5kg



Compact Wheelbase
Rear wheels are positioned for maximum clearance during entry and exit from the frame



Easy Folding
The walker folds easily for convenient transport and compact storage



Removable Storage
Zip bag secures up to 5kg of essentials

Great for Outdoors The large 10" wheels provide stability and smooth movement across a variety of terrains

Soft-touch Brakes Lock the rear wheels to enhance safety and stability during use

Adjustable Handles Allows for personalised positioning, helping to maintain proper posture

Ergonomic Handles These offer a comfortable grip during extended use. The flat-top design provides a stable resting place for arms while seated

Bonus Walking cane holder

SWL
150kg

SWL
23.6st

1 YEAR WARRANTY


Black

	Overall Length	Overall Width	Overall Height	Seat Width	Seat Height	Folded Width	Wheel Size Front & Rear	Minimum Functional Weight	Typical Functional Weight**
	mm / "	mm / "	mm / "	mm / "	mm / "	mm / "	mm / "	mm / "	kg / lbs
Black - Medium WAF705380	710	625	860 - 945	460	540	245	250 & 200	5.4	6.35
	28	24.5	34 - 37	18	21	9.6	10 & 8	12	14
Black - Tall WAF705390	710	625	860 - 1025	460	620	245	250 & 200	5.5	6.4
	28	24.5	34 - 37	18	24.4	9.6	10 & 8	12	14

VOGUE LIGHTWEIGHT 2 SEAT WALKER

Lightweight Aluminium Body Improves manoeuvrability and portability

This modern European designed seat walker has a clean and crisp frame design and is available in a selection of eye-catching colours and prints to really stand-out while keeping you on the move.

Just Under
6kg



Compact Turning Circle
For improved mobility in tight spaces



Ergonomic Handles
With personalised height adjustment to promote comfortable and correct posture



Handbrakes Lock
For user safety and stability


Lightweight Aluminium Body Improves manoeuvrability and portability

Simple Folding Mechanism Locks for easy transportation and storage

Removable Storage Zip bag secures up to 5kg of essentials

Deep Seat and Wide Backrest For improved comfort and support

Bonus Walking cane holder



Leopard
WAF705310LE



Flower
WAF705310FL



Zebra
WAF705310ZB



Red
WAF705310RE



Champagne
WAF705310CH

SWL
150kg

SWL
23.6st

1 YEAR WARRANTY

Overall Length mm / "	Overall Width mm / "	Overall Height mm / "	Seat Width mm / "	Seat Height mm / "	Folded Width mm / "	Wheel Size mm / "	Minimum Functional Weight kg / lbs	Typical Functional Weight** kg / lbs
655	620	780 - 920	465	540	220	200	5.4	6.45
26	24.5	31 - 36	18	21	8.6	8	11.7	14.2

VOGUE ADVENTURE SEAT WALKER

Equipped with Extra-large All-terrain Wheels Built for travel beyond the everyday

Rugged and durable, the Adventure Seat Walker is the perfect all terrain walker. Super-sized front wheels can manage any surface, perfect for the adventurous user.



Large Multi-terrain Wheels
For durability and stability on and off road



Simple Folding Mechanism
Locks for easy transportation and storage



Removable Storage
Zip bag secures up to 5kg of essentials

Ergonomic Design Offers personalised height adjustment which promotes correct posture

Heavy-duty Frame and Wider Seat Provide enhanced support for users up to 200kg

Handbrakes Lock For user safety and stability

Compact Turning Circle for improved mobility in tight spaces

Bonus Walking cane holder

SWL
200kg

SWL
31.5st

1 YEAR WARRANTY

Overall Length mm / "	Overall Width mm / "	Overall Height mm / "	Seat Width mm / "	Seat Height mm / "	Folded Width mm / "	Wheel Size Front & Rear mm / "	Minimum Functional Weight kg / lbs	Typical Functional Weight** kg / lbs
770	695	900 - 1035	540	530	245	279 & 225	8.8	10
30.3	27.4	35.4 - 40.7	21.3	20.9	9.6	11 & 8.8	19.4	22




Black
WAF705450

VOGUE FOREARM SEAT WALKER

Lightweight Aluminium Body Improves manoeuvrability and portability

For those who need a little more support, the Vogue Forearm Seat Walker is for you. Soft and contoured forearm pads ensure maximum comfort and a smooth ride.



Smart Posture System
With adjustable forearm supports across height, depth and angle

Simple Folding Mechanism
Folds to 225mm

Removable Storage
Zip bag secures up to 5kg of essentials

Luxury Look and Feel Danish design, profiled tubing, moulded fittings, ergonomic handles, and refined colours

Lightweight and Ergonomic With small turning circle for excellent manoeuvrability

AutoLocate Smooth Height Adjustment That automatically locates to pre-set heights, making adjustment simple and accurate

Locks when Folded Makes transportation and storage easier

Advanced Locking Handbrakes For increased user safety when walking and sitting

Unbreakable TPE Super-soft Wheels Offers high durability and a smooth ride for increased comfort

SWL
150kg

SWL
23.6st

1 YEAR WARRANTY

White / Black
WAF705400

Overall Length mm / "	Overall Width mm / "	Overall Height mm / "	Seat Width mm / "	Seat Height mm / "	Folded Width mm / "	Wheel Size Front & Rear mm / "	Minimum Functional Weight kg / lbs	Typical Functional Weight** kg / lbs
655	620	900 - 1170	441	620	225	200	9.2	10
26	24	40 - 50	17	24.5	9	8	20	22

VOGUE INDOOR WALKER

Lightweight Aluminium Body For easy mobility and meal transport around the home

Ideal for use within the home, the Indoor Walker is perfect for portable dining and room to room mobility support.



Simple Folding Mechanism
Folds to 270mm, the slimline design reduces the footprint when folded

No Lifting and Spillage
Meals can be collected and brought to the desired location without the need for lifting or carrying. The tray is also removable

Removable Storage Bag
Utensils and other items can be stowed securely away, ideal for transporting personal belongings around the house

Luxury Look and Feel Danish design, profiled tubing, moulded fittings, ergonomic handles, and refined colours

Lightweight and Ergonomic Design With small turning circle for excellent manoeuvrability

Sitting Down for a Meal When the user is sitting, the walker can be turned around and positioned comfortably in front of them, to enjoy their meal

Brakes The ergonomic 'Soft Touch' lockable design allows the walker to be secured in place with minimal hand pressure

SWL
110kg

SWL
17.3st

1 YEAR WARRANTY

Red / Black
WAF705100RE

Champagne / Black
WAF705100CH

Overall Length mm / "	Overall Width mm / "	Overall Height mm / "	Folded Width mm / "	Wheel Size mm / "	Minimum Functional Weight kg / lbs	Typical Functional Weight** kg / lbs	Tray SWL kg / lbs
650	540	860 - 970	270	152	4.4	5.1	5
25.6	21.3	33.9 - 38.2	10.7	6	9.7	11.2	11



VOGUE WALKER ACCESSORIES



Oxygen Bottle Holder Bag*
WAA695600



LED Night Light
WAA695630



Height Adjustable IV Pole
WAA695651



Single Hand Brake System*
WAA695640



Cup Holder**
SCA669200



Phone Holder**
SCA669210

*Single Hand Brake System and Oxygen Bottle Holder Bag are not compatible with the Vogue Indoor Walker **Non-Aspire branded products

Mobility Accessory Bags

Keep personal items close and secure when you're on the move. Ideal for storing your everyday essentials, including your phone, wallet, or water bottle with a 5kg bag capacity.



Leopard WAA696500



Floral WAA696510



Plaid WAA696520



Blue WAA696530

