



STAY IN THE HOME  
**YOU LOVE**



## BRUNO'S BEST-SELLING STAIRLIFT, THE ELAN

Bruno's most popular stairlift can help restore full access to your home. Designed for straight indoor staircases, the Elan offers a mix of high-quality standard functionality, unsurpassed reliability, and ease-of-use design.

With high quality and dependability, Bruno helps you stay in your home with comfort and ease.



Verified Customer

"Since we got the stairlift, we don't know what we did without it for so long. We are very pleased with it."

Joseph L.





## DEPENDABLE QUALITY. QUIET RIDE.

**Easy Operation:** Push the armrest paddle control to travel up and down the stairs. Two wireless remote controls also included.

**Space Saving:** Narrow vertical rail installs close to the wall to preserve maximum stairway space.

**Sleek & Clean:** The Elan's vertical rail with hidden gear rack gives a clean appearance and reduces exposure to dirt and grease.

**Dependable:** Battery-powered operation ensures performance even in a power outage. Stairlift continually charges along entire rail.

**Compact:** Fold the Elan's arms, seat and footrest to save space when stairlift is not in use.

**Safe:** Seat swivels up to 90° away from steps, then locks in place for safe exit. Seat belt and obstruction sensors ensure safety.

[Learn More](#)  
[bruno.com/aus](https://bruno.com/aus)



Verified Customer

"I love the Bruno stairlift! It takes up much less space than I expected."


Mary Jo C.



## AVAILABLE OPTIONS

Add a folding rail for a narrow hallway application or when a doorway is at the bottom of the stairs. Manual operation or push-button automatic.

Convert the manual Elan swivel seat to optional power swivel and the chair will automatically rotate with the push of the armrest control. Further simplify your Elan with a power folding footrest. Fold the seat and the footrest automatically folds, too. Great for people with back issues!



Verified Customer

"Excellent quality, runs smoothly,  
money well spent."

Priscilla M.



## TECH SPECS & WARRANTY

**Lift Capacity:**

136 kg (300 lb)

**Seat Pad Size:**

400 mm wide x 426 mm deep (15.75 in x 16.75)

**Width Between Arms:**

514 mm (20.25 in)

**Footrest Size:**

267 mm long x 279 mm wide (10.5 in x 11 in)

**Vertical Rail:**

Installs within 127 mm (5 in) of the wall



## ELAN WARRANTY

**Lifetime Warranty on Major Components** including motor, gearbox, drivetrain and stationary rail system, for the lifetime of the Bruno Stairlift as long as the original purchaser owns the stairlift.

**2 Year Warranty** on all components after date of installation by an Authorised Bruno Dealer.

**30 Days Labour Coverage** after date of installation by an Authorised Bruno Dealer for any manufacturing defect.

## ABOUT BRUNO

**Veteran-Founded | Family-Owned**

Bruno makes stairlifts and vertical platform lifts to help people better access homes, public buildings and businesses, and scooter lifts to assist people in easily transporting their mobility devices.

