



ACCESS ALL LEVELS OF
YOUR HOME



ENJOY BRUNO'S SIGNATURE STAIRLIFT

Experience the ultra-smooth, quiet ride of a Bruno Elite curved stairlift. Bruno is the only manufacturer who handcrafts the rail to ensure a precision fit that artfully flows into every corner and maximises open space on the steps.

With heavy-duty capacity, high performance engineering and a Limited Lifetime Warranty, you gain the peace of mind that a Bruno Elite curved stairlift will allow you to stay safely in your home long into the future.



3 TOP CONSIDERATIONS

- ✓ Handcrafted by skilled welders
- ✓ 181kg (400 lb) Lift Capacity
- ✓ Limited Lifetime Warranty





CUSTOM FIT. HIGH-PERFORMANCE.

Handcrafted: Elegant, handmade rail delivers a precision fit and smooth ride while the covered gear rack gives clean, graceful appearance.

Easy to Use: Armrest control glides stairlift up and down. Two wireless call/send remotes are also included.

Adjustable: Change seat and footrest height and space between arms for optimal fit for your body.

Compact: Vertical, custom rail leaves maximum open space on steps. Fold the arms, seat and footrest when not in use to save even more space.

Safe: Seat swivels 90° at top of steps and locks in place for safe exit away from stairs. Seat belt and obstruction sensors further ensure safety.

Dependable: Battery-powered operation ensures performance, even in the event of a power outage.

Learn More
bruno.com/aus

AVAILABLE OPTIONS

Add extended rail positions to park the stairlift safely away from the steps or out of view and an optional mid-park charging station for stairs with middle landings.

Get an extra-wide seat or larger footrest to comfortably accommodate larger stairlift riders.

Eliminate the need to lean over with a power folding footrest. Fold the seat and the footrest will automatically fold, too. Convert from manual to a power swivel seat and the chair will automatically rotate with the push of a button.

Select from four additional seat fabric colour options to blend your stairlift perfectly into your home's decor.



Not all options are shown.



Verified Customer

"I love my Bruno stairlift. It works so smoothly... It is truly a blessing."

Sharon W.





TECH SPECS

Lift Capacity:

181 kg (400 lb)

Standard Seat Size:

457 mm wide x 381 mm deep (18 in x 15 in)

Larger Seat Option:

572 mm wide x 419 mm deep (22.5 in x 16.5 in)

Standard Footrest Size:

254 mm long x 229 mm wide (10 in x 9 in)

Larger Footrest Option:

305 mm long x 279 mm wide (12 in x 11 in)

Curved Rail:

Installs within 190 mm (7.5 in) of the wall

WARRANTY

Lifetime Warranty on Major Components for the lifetime of the Bruno Stairlift as long as the original purchaser owns the stairlift.

2 Year Warranty on all components after date of installation by an Authorised Bruno Dealer.

30 Days Labour Coverage after date of installation by an Authorised Bruno Dealer for any manufacturing defect.

ABOUT BRUNO

Veteran-Founded | Family-Owned

Bruno makes stairlifts and vertical platform lifts to help people better access homes, public buildings and businesses, and scooter lifts to assist people in easily transporting their mobility devices.

