

# ICON 30 FAF

FOR OUTDOOR ACTIVITIES



REHASENSE



Exclusively distributed by

**aidacare™**  
HEALTHCARE EQUIPMENT

# ICON 30 FAF

## THE LIGHTWEIGHT WHEELCHAIR

ICON 30 FAF is a stable, foldable wheelchair specifically designed for outdoor activities. It's fixed anterior frame (FAF) and multiple adjustment options allow it to be tailored specifically for the users needs to maximise performance and comfort. Constructed from high-strength aluminium tubing with reinforced side frames, it provides a stable and safe connection for PAWS and Eco Assist Power Add-ons.



**Foldable Rear Handles**  
Enables the chair to be propelled by a carer or family member



**Camber**  
Up to 6° of rear wheel camber



**Drop-in Armrests**  
Facilitate side transfers



**Footplates**  
Height adjustable to suit multiple user heights

### REAR WHEEL OPTIONS

Sold separately



**MWT481800**  
FATBOY 20" Smooth Tread



**MWT481810**  
FATBOY 20" Offroad Tread



**MWT481820**  
EASYSRIDER 24" Smooth Tread



**MWT481830**  
FATMAN 24" Smooth Tread



**2 YEAR WARRANTY**

ICON 30 FAF	Seat Width mm	Seat Depth mm	Anterior Seat Height mm	Posterior Seat Height mm	Overall Width mm	Backrest Height mm	Rear Wheel Size "	Front Castor Size "	Min Chair Weight kg	SWL kg
MWC451400	400	400	435 - 485	390 - 455	630	310 - 340	24	5	from 13.6	125
MWC451410	400	450	435 - 485	390 - 455	630	310 - 340	24	5	from 13.6	125
MWC451420	450	400	435 - 485	390 - 455	680	310 - 340	24	5	from 13.6	125
MWC451430	450	450	435 - 485	390 - 455	680	310 - 340	24	5	from 13.6	125