

ICON 60

THE ACTIVE WHEELCHAIR



REHASENSE



Exclusively distributed by

aidacare™
HEALTHCARE EQUIPMENT

ICON 60

THE ACTIVE WHEELCHAIR

ICON 60 is an adjustable, lightweight wheelchair made from aluminium. To make the chair more sporty or comfortable, users can make many adjustments including the seat depth, front and rear seat height, seat angle, and seat position to suit their needs and preferences. ICON 60 also features a backrest that can tilt forward up to 20°, or lean back up to 10°, which can be folded down completely to transport easily.



ICON 60 with Eco Assist
Sold separately: WCA872000

- Active wheelchair in ‘open frame’ design
- Multi-adjustable: Seat depth, front and rear seat heights, and centre of gravity
- Angle adjustable front forks



ICON 60 with PAWS City 14" Auto-clamp
Standard Control Sold separately: WCA872112

- Angle, height, and tension adjustable backrest
- Breathable 3D mesh fabric
- Open footrest design

REAR WHEEL OPTIONS

Sold separately



MWT481800 FATBOY
20" Smooth Tread



MWT481810 FATBOY
20" Offroad Tread



MWT481820 EASYRIDER
24" Smooth Tread



MWT481830 FATMAN
24" Smooth Tread



Closed backrest with foldable handles for attendant propelling support



Carbon fibre side panels

2 YEAR WARRANTY

ICON 60	Seat Width mm	Seat Depth mm	Anterior Seat Height mm	Posterior Seat Height mm	Backrest Height mm	Rear Wheel Size "	Front Castor Size "	Min Chair Weight kg	SWL kg
MWC451500	360	420 - 465	490 - 510	490 - 510	300	24	5	From 6.9	125
MWC451510	380	420 - 465	490 - 510	490 - 510	300	24	5	From 6.9	125
MWC451520	400	420 - 465	490 - 510	490 - 510	300	24	5	From 6.9	125
MWC451530	420	420 - 465	490 - 510	490 - 510	300	24	5	From 6.9	125
MWC451540	440	420 - 465	490 - 510	490 - 510	300	24	5	From 6.9	125
MWC451550	460	420 - 465	490 - 510	490 - 510	300	24	5	From 6.9	125
MWC451560	480	420 - 465	490 - 510	490 - 510	300	24	5	From 6.9	125