

## Aidacare FL250 Floorline Bed



- Under Bed Clearance - While the FL250 achieves a low height of 110 mm, full clearance under the bed is still maintained, minimising the risk of foot or cord impingement
- SmartDescent™ Automatically slows descent speed over the final 300 mm of lowering to the floor level, allowing possible impingements to be identified and moved
- Hidden Bed Mechanisms: The bed mechanisms are attached under the platform and are not visible. This helps the room feel less clinical
- Low Height 11 cm: Adding a 15 cm mattress gives you a total height below 30 cm thus significantly reducing out-of-bed fall injury risks
- Auto-Regression on back & knee rests: Reduces skin shear forces during operation thus improving clinical outcomes and reducing the need for expensive mattress solutions.
- High SWL 250 kg & base-extension to 100 & 200 mm: The same bed platform can meet a wider profile of residents
- BrakeFirst™ Automatic Safety Brake Activation at lowest height: Reduces the risk of injury during bed exit in the morning as the bed will lock when at its lowest height
- Single Button Seat Position: Easy transfer of weight from lower back to bottom for pressure-care and seated positioning also helps improves lung function for resident
- ExitSafe™ Single Button Ingress/Exit height: Your staff don't have to think about what height to set the bed. Via one button the vast majority of needs are met
- Attendant Control Panel - Allows the caregiver to lock out certain functions and also ensures bed operation is not just restricted to the handset
- Grey Colour Scheme & Narrow-profile frame: Keeps the focus on bedding and furnishings/finishings rather than the bed-frame

✓ Contracted DVA product

✓ Contract approved

### Support @ Home

**Function:** Beds and bed equipment

**Category:** Prescribed

Click [here](#) for more information

### Compliance

IEC 60601-1:2005

IEC 60601-2-52:2005

### Accessories

Click on the cart icon to display the range of optional accessories

- Underbed light standard: Helps to reduce incidence of resident and caregiver falls
- Rubber & Nylon castors: Castors remain on the ground during braking, thus not damaging floors
- King Single Extension Kit (optional retrofit): Allows you to change the profile of your beds (with minimal cost) to meet resident profile changes
- KTL coated anti-corrosion steel: Improved durability for lower lifecycle cost
- Low-maintenance Braking System: Fewer headaches for your MO's!

## Warranty Information

- 10 Years Frame
- 3 Years Electrics
- 1 Year Handset

| Code      | Display Name  | Total Length | Total Width | Total Height | Safe Working Load | Colour        |
|-----------|---|--------------|-------------|--------------|-------------------|---------------|
| BEB046100 | Aidacare FL250 Floorline Bed - Single (Bed ends sep) - Incl: UBL & Extension      | 2000 mm      | 900 mm      | 110 – 740 mm | 250 kg            | Sublime Teak  |
| BEB046150 | Aidacare FL250 Floorline Bed - King Single (Bed ends sep) - Incl: UBL & Extension | 2000 mm      | 1070 mm     | 110 - 740 mm | 250 kg            | Milano Walnut |