

Aspire Transfer Bench with Rotating Seat



- Comfortable padded backrest, seat and arms offer a warmer sitting surface
- Height adjustable legs with 12.5 mm increments suit a range of bath to floor heights
- Contoured seat with a safety belt, slides smoothly and is easy to clean
- 360° rotating seat (with lockable positions every 90°) suits left and right bath orientation
- Slip-resistant suction rubber tips for added stability
- Lightweight aluminium construction for ease of use and portability
- SWL - 136 kg

✓ Contract approved

Support @ Home

Function: Washing, bathing and showering

Category: Prescribed

Click [here](#) for more information

Warranty Information

- 1 Year Warranty



Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load	Material	Unit Weight	Seat Depth	Seat Width	Seat To Floor Height	Seat To Footrest
BTS150310	Aspire Transfer Bench with Rotating Seat	475 mm	950 mm	795 - 905mm	136 kg	Aluminium Frame	8.1 kg	345 mm	425 mm	505 - 615 mm	505 - 615 mm