

Aspire PosturePro Chair



Optimised Seating for Comfort and Support -The PosturePro is designed for users requiring medium to high postural support and additional comfort during prolonged seating. Combining memory foam, air bladders, and multi-stretch fabric, it ensures pressure relief while maintaining user comfort.

- **Customisable Backrest - Integrated Laterals:** Offers side support to promote upright posture and prevent lateral leaning.
- **Adjustable Pocketed Waterfall Design:** Enables personalised back support by adjusting to the user's specific back shape and requirements.
- **Tool-Free Seating Adjustment - Seat Depth and Width Configuration:** Allows therapists or carers to adjust seat dimensions quickly and without tools, ensuring the user achieves a comfortable and appropriate sitting position.
- **Adaptable Footrest - Depth and Angle Adjustability:** Customises footrest depth and angle to meet the user's positioning and weight-bearing needs, aiding sit-to-stand transfers.
- **Stand Air Compatibility - Safe Standing Transfers:** Integrates with the Aspire Star Mover to facilitate safe, efficient standing transfers, enhancing user mobility and independence.
- **Arm Height Adjustability - Height and Depth Adjustments:** Adjusts armrests in both height and depth to ensure proper arm loading, providing additional comfort and support.
- **Tilt-in-Space - Infinite Positioning Options:** Enables tilting and reclining at various angles, making repositioning easier and enhancing comfort.
- [Watch instructional adjustment video.](#)
- Frame - 2 Years
- Components - 1 Year
- Fabric - 1 Year

✓ Contract approved

Support @ Home

Function: Sitting furniture

Category: Under advice

Click [here](#) for more information

Code	Display Name	Total Width	Total Height	Safe Working Load	Seat Depth	Seat Width	Seat To Floor Height	Seat To Footrest	Backrest Height
CHP231600	Aspire PosturePro Chair	650 - 850 mm	1300 - 1360 mm	140 kg	450 - 540 mm	330 - 530 mm	530 - 590 mm	530 - 590 mm	770 mm