

## Telescopic IV Pole Malvestio Botero Bariatric Ward Bed



The specific bariatric design focus offers width and length adjustable deck to manage changing needs in both ward and respite environments.

- Low height of 380 mm
- Antibacterial powder coating for improved infection control management
- Backrest translation (autoregression) reduces compression of the thoracic and abdominal areas as the backrest is raised and alleviates shear and pressure for improved skin integrity
- Compatible with all leading bed movers
- Attendant control panel with intuitive pre-programmed single touch functions for speed and simplicity
- 150 mm TENTE Castors
- Central locking castors with integrated directional lock
- Full Trendelenburg functionality
- Integrated folding side rails
- Removable width adjustable head and foot panels
- Back rest angle 0 to 70 degrees
- Knee bend angle 0 to 30 degrees
- Trendelenburg -17 to 17 degrees

### Warranty Information

- 2 Years



Contract approved

### Support @ Home

**Function:** Beds and bed equipment

**Category:** Prescribed

Click [here](#) for more information

### Compliance

IEC 60601-2-52:2009

Code	Display Name
BEA825101	Telescopic IV Pole - Malvestio Botero Bariatric Ward Bed