


## Aspire ComfiMotion Allure Footboard



Aspire ComfiMotion headboards and footboards blend style with function, helping to create a dignified and practical sleeping environment. Whether freestanding or fixed, each variant addresses specific clinical and aesthetic needs—improving comfort for users and making life easier for carers.



The footboard completes the bed in both style and function, compatible with **ComfiMotion ActiveCare beds only.**

- Helps complete the bed end and keep bedding in place during profiling, reducing the need for carers to continually re-adjust bedding
- Provides a practical accessory mounting surface (e.g. mattress pump)
- Portobello grey upholstery to match bed and headboards
- Secure and sturdy bolt-on attachment to bed frame
- Widths to match bed sizes, all 610mm high and 60mm deep

### Warranty Information

- 1 Year

Code	Display Name	Total Length	Total Width	Total Height	Colour
BEA830166	Aspire ComfiMotion Footboard - Allure - Portobello Grey - Queen	60 mm	1520 mm	610 mm	Portobello Grey
BEA830164	Aspire ComfiMotion Footboard - Allure - Portobello Grey - Long Double	60 mm	1370 mm	610 mm	Portobello Grey
BEA830162	Aspire ComfiMotion Footboard - Allure - Portobello Grey - King Single	60 mm	1070 mm	610 mm	Portobello Grey
BEA830160	Aspire ComfiMotion Footboard - Allure - Portobello Grey - Long Single	60 mm	920 mm	610 mm	Portobello Grey

