

Malvestio Delta 4 Ward Bed



The DELTA4 Ward Bed is the new generation of beds developed according to the needs of healthcare professionals. The DELTA4 bed becomes an active platform that simplifies and maximises the performances of healthcare professionals.

- **INCENTRO System** - the combined back and pelvic translation nullifies horizontal forces, prevents slipping and keeps the body at a 30° angle, optimal for ventilation
- **Safe Exit Lights** - 2 x LED lights signal whether the bed is at a safe exit height, promoting safe egress for the patient (orange is unsafe and green is safe)
- **Safe Exit control** moves the platform to the minimum height and lifts the backrest up to facilitate the patient to get up from the bed
- **Auto Regression** - backrest automatically shifts backwards thus increasing the space on the pelvic area by creating an extra space of 13 cm
- **Biocote** - metal components are epoxy powder coated and finished with an antibacterial additive for improved infection control
- **DeltaDrive** - electric centre 5th wheel with integrated controls in headboard provides increased manoeuvrability for hospital staff (available on BEB050150 only)
- **Weigh Scale** - integrated system allows for precision weight measurements that can be set to automatically weigh a patient every 24 hours (available on BEB050150 only)
- **4 sector safety sides** with integrated electric controls which allows full autonomy of movement to the patient
- **CPR control** on both sides allows the bed platform to be brought to the horizontal position immediately accessible
- **Mattress platform** can be easily removed, all sections are designed for easy cleaning
- **Automatic self-brake** once connected to the power net a sound alarm alerts

✓ Contract approved

Support @ Home

Function: Beds and bed equipment

Category: Prescribed

Click [here](#) for more information

the nursing staff about the unbraked bed and after 20 seconds the brake is automatically activated

- **X-ray tray** together with HPL backrest allows extremely easy x-ray operations to be carried out without moving the patient
- **StaminaLift Compatibility** - allows for use with existing hospital systems (available upon request for BEB050140 only)

Warranty Information

- 10 Years Frame
- 2 Years Electrics
- 1 Year Handset

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load	Unit Weight
BEB050140	Delta 4 Ward Bed	2200 mm	990 mm	390 - 790 mm	300 kg	125 kg