

## Aspire Signature 2 CUSTOM Lift Recline Chair - Dual Action



Customise your Aspire Signature 2 - Dual Action Lift Recline Chair with upholstery to match your style and needs.



- **Healthcare Vinyl** The Lustrell Ranger range is high performance vinyl with HealthGuard Antibacterial and Antimicrobial treatment to assist with infection control and hygiene Commercial quality fabric for long lasting durability
- **Microsuede** Suede like texture combined with high performance characteristics ensures the Eastwood fabric range remains easy to clean and maintain
- **Leather** La Casa is a luxurious artisan leather with a soft glove-like texture, will retain their unique suppleness, provided they are regularly cared for
- **Soft Touch Fabric** The soft touch weave and diverse colour palette makes the Mazza and Manisa ranges a highly practical and functional solution
- Commercial quality fabric for long lasting durability
- Custom order lead time varies, approximately 5 - 6 weeks depending on fabric availability
- Custom orders are subject to terms and conditions please see page two on the order form below
- All colour swatches shown may not be a true representation of the actual product and should be used as a guide only
- **Raphael Powerlift Recline Chair Features & Benefits:**
  - **Dual Action:** Independent control of legrest and recline functions
  - **Multiple Sizes:** Available in small, medium, and large sizes to ensure perfect user fit
  - **Leg Extension:** Helps to comfortably support lower legs for a wide range of user heights..
  - **Waterfall Backrest:** Pocketed design allows filling to be adjusted to

suit unique user postural requirements

- **All-Day Battery Backup:** Ensures safety and user confidence during power outages
- **Easy to Move:** Four swivel castors facilitate movement when the chair needs to be moved
- **Depth Adjustable Seat:** A choice of two depth settings in 1in (25mm) increments for custom comfort and positioning
- **Integrated USB Charging:** Conveniently located at the end of the handset, which is backlit

## Warranty Information

- 2 Years

Code	Display Name	Total Width	Total Height	Safe Working Load	Material	Seat Depth	Seat Width	Seat To Floor Height	Seat To Footrest	Backrest Height
CHP228280VL	Signature 2 CUSTOM Lift Recline Chair - Dual Action - Large - Healthcare Vinyl (choose colour)	820mm	1110mm	155kg	Vinyl	510-535mm	540mm	500mm	500mm	650mm
CHP228280MS	Signature 2 CUSTOM Lift Recline Chair - Dual Action - Large - Microsuede (choose colour)	820mm	1110mm	155kg	Microsuede	510-535mm	540mm	500mm	500mm	650mm
CHP228280LT	Signature 2 CUSTOM Lift Recline Chair - Dual Action - Large - Leather (choose colour)	820mm	1110mm	155kg	Leather	510-535mm	540mm	500mm	500mm	650mm

Code	Display Name	Total Width	Total Height	Safe Working Load	Material	Seat Depth	Seat Width	Seat To Floor Height	Seat To Footrest	Backrest Height
CHP228280FA	Signature 2 CUSTOM Lift Recline Chair - Dual Action - Large - Soft Touch Fabric (choose colour)	820mm	1110mm	155kg	Soft Touch Fabric	510-535mm	540mm	500mm	500mm	650mm
CHP228250VL	Signature 2 CUSTOM Lift Recline Chair - Dual Action - Medium - Healthcare Vinyl (choose colour)	780mm	1050mm	155kg	Vinyl	490-515mm	510mm	480mm	480mm	620mm
CHP228250MS	Signature 2 CUSTOM Lift Recline Chair - Dual Action - Medium - Microsuede (choose colour)	780mm	1050mm	155kg	Microsuede	490-515mm	510mm	480mm	480mm	620mm
CHP228250LT	Signature 2 CUSTOM Lift Recline Chair - Dual Action - Medium - Leather (choose colour)	780mm	1050mm	155kg	Leather	490-515mm	510mm	480mm	480mm	620mm
CHP228250FA	Signature 2 CUSTOM Lift Recline Chair - Dual Action - Medium - Soft Touch Fabric (choose colour)	780mm	1050mm	155kg	Soft Touch Fabric	490-515mm	510mm	480mm	480mm	620mm
CHP228220VL	Signature 2 CUSTOM Lift Recline Chair - Dual Action - Small - Healthcare Vinyl (choose colour)	740mm	980mm	155kg	Vinyl	460-385mm	480mm	480mm	480mm	580mm
CHP228220MS	Signature 2 CUSTOM Lift Recline Chair - Dual Action - Small - Microsuede (choose colour)	740mm	980mm	155kg	Microsuede	460-385mm	480mm	480mm	480mm	580mm
CHP228220LT	Signature 2 CUSTOM Lift Recline Chair - Dual Action - Small - Leather (choose colour)	740mm	980mm	155kg	Leather	460-385mm	480mm	480mm	480mm	580mm

Code	Display Name	Total Width	Total Height	Safe Working Load	Material	Seat Depth	Seat Width	Seat To Floor Height	Seat To Footrest	Backrest Height
	colour)									
CHP228220FA	Signature 2 CUSTOM Lift Recline Chair - Dual Action - Small - Soft Touch Fabric (choose colour)	740mm	980mm	155kg	Soft Touch Fabric	460-385mm	480mm	480mm	480mm	580mm