

## Malvestio Botero Static Mattress



- Two layer foam design with high density base and castellated surface caters to bariatric patients
- Multi-stretch cover with welded seams for stringent infection control
- Strong high density side walls improve patient transfer efficiency while cut back hinge points allow the mattress to easily move with the mattress platform

### Warranty Information

- Foam: 5 Years
- Cover: 2 Years



Contract approved

### Support @ Home

**Function:** Manage tissue integrity

**Category:** Prescribed

Click [here](#) for more information

### Compliance

BS7177:2008

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load
BEM050200	Botero Static Mattress	2000 mm	1220 mm	200 mm	400 kg