

## Malvestio VIVO ICU Bed



The VIVO bed was designed with patient and care staff welfare in mind. It offers integrated features to ease care staff burden whilst creating a comfortable and effective patient support surface.

- Lateral tilting patient surface
- Alternating air mattress with full bed integration to maximise therapeutic benefit
- Patient Anti Slipping System (PAS) offsets the patient movement towards the foot end of the bed during backrest adjustment
- Radio transparent patient surface for fluoroscopic examinations
- Touch screen panel in the foot end allows for the bed and integrated air mattress functions to be controlled with ease
- Integrated bed extension to suit a range of patient heights
- Intuitive control panels within the gas assisted folding side rails
- Safety alarm system ensures the base of the bed is clear before any movements, preventing entrapment and damage
- Backrest angle of 0 to 68 degrees
- Knee bend Angle of 0 to 25 degrees
- Trendelenburg - 15 to 15 degrees

### Warranty Information

- 2 Years



Contract approved

### Support @ Home

**Function:** Beds and bed equipment

**Category:** Prescribed

Click [here](#) for more information

### Compliance

IEC 60601-2-52:2009

| Code      | Display Name            | Total Length      | Total Width | Total Height | Safe Working Load |
|-----------|-------------------------|-------------------|-------------|--------------|-------------------|
| BEB050100 | Malvestion VIVO ICU Bed | 2190 mm - 2420 mm | 1080 mm     | 480 - 820 mm | 260 kg            |